

# FLEX BENCH

## SERIES

**FLEX BENCH SERIES** Our highly adaptable flex benches are offered in many configurations and suitable for a wide variety of laboratory functions. Benches are equipped with heavy-duty casters and designed to accommodate a full range of overhead and suspended storage components. The integrated supply carrier provides easy access to electric and data services at the benchtop.



Height Adjustable 30-36"

Heavy-Duty Steel Frame

Fixed or Mobile Base Style



Easy Access to Services

# FLEXBENCH

## FEATURES

- + **Height Adjustable** [30-36"// sit to stand table position]
- + **Fixed Style or Mobile Base** [4" heavy duty locking casters]
- + **Powder Coated Steel Frame** [available in a variety of colors]
- + **Durable Work Surface Options** [phenolic resin or plastic laminate]
- + **Service Supply Modules** [easy access to H/C water, gas, electricity & data]
- + **Under-Bench Storage** [optional drawer, cupboard or combination storage]
- + **Adjustable Reagent Shelving** [with optional LED task lighting]



## TABLE FRAMES



### STEEL COLORS [60% Gloss]

WHITE | BLACK | STORM GRAY | BEIGE | BLUE | DARK BLUE



## WORK SURFACES



### PHENOLIC RESIN COLORS

TITAN WHITE | BLACK



### PLASTIC LAMINATE COLORS

WHITE | BLACK | CLOUDY GRAY | BEIGE | STORM GRAY

## DIMENSIONS

ITEM #	TABLETOP	HEIGHT ADJUST	FRAME	WORK SURFACE OPTIONS
LS-KY2424	24"D X 24"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2430	24"D X 30"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2436	24"D X 36"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2448	24"D X 48"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2460	24"D X 60"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2472	24"D X 72"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2496	24"D X 96"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate

▶ **ALSO AVAILABLE** 30" | 36" | 48"D