

INSTRUCTOR DESK

DEMONSTRATION CENTER

INSTRUCTOR DESKS accommodate storage needs and demonstration requirements in one unit. They are designed to serve as a teaching aid for lab experiments and a standard workstation during lecture time. Optional features include a height adjustable worktop, mobile base, integrated computer technology and services for gas, electric and plumbing. Instructor Desks are furnished with one upright rod assembly and can be configured with a sink on either left or right side.

FEATURES

Multipurpose Design

Height Adjustable Feature
available for some units

Chem-Resistant Epoxy Resin Top

Top Sizes from 5' to 9' in Length

Integrated Computer Technology

Integrated Service Options



STANDARD DEMONSTRATION UNIT



PORTABLE DEMONSTRATION UNITS

LAB TABLES

SOLID WOOD FRAME CONSTRUCTION [CHEMICAL RESISTANT - UV FINISH]

5+ WOOD FINISH OPTIONS

FIXED STYLE OR MOBILE BASE

DURABLE WORK SURFACE OPTIONS [EPOXY OR PHENOLIC RESIN]

METAL WIRE PULLS STANDARD

PROTECTIVE RUBBER BOOTS

HEAVY-DUTY LOCKABLE CASTERS



TABLE FRAMES



STEEL COLORS [60% Gloss]

WHITE | BLACK | STORM GRAY | BEIGE | BLUE | DARK BLUE



WORK SURFACES



PHENOLIC RESIN COLORS

TITAN WHITE | BLACK



PLASTIC LAMINATE COLORS

WHITE | BLACK | CLOUDY GRAY | BEIGE | STORM GRAY

DIMENSIONS

ITEM #	TABLETOP	HEIGHT ADJUST	FRAME	WORK SURFACE OPTIONS
LS-KY2424	24"D X 24"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2430	24"D X 30"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2436	24"D X 36"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2448	24"D X 48"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2460	24"D X 60"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2472	24"D X 72"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2496	24"D X 96"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate

▶ ALSO AVAILABLE 30" | 36" | 48"D