


LAB TABLES

// SYSTEMS / COMPONENTS / ACCESSORIES




FEATURES

METAL FRAME

- + 16" or 18" gauge steel frame construction
- + Work surface in epoxy or phenolic resin
- + Powder coated steel finish
- + Casters / adjustable height legs optional
- + Metal wire pulls standard

COLORS



PETAL WHITE | LIGHT ALMOND | CHAMPAGNE | OFF WHITE | MODEL GRAY
BLACK | BANKER'S GRAY | PEWTER | HUNTER GREEN | STONE GRAY
FAWN BEIGE | MEDIUM BLUE | WEDGEWOOD BLUE | BURGUNDY

WOOD FRAME

- + Solid wood frame construction
- + Chemical resistant, UV finish
- + Work surface in epoxy or phenolic resin
- + Available in 5+ finishes
- + Metal wire pulls standard


FINISHES




NATURAL OAK | CINNAMON TOAST | FALL HARVEST | GINGER | NORTHERN OAK

STANDARD DIMENSIONS

D 24" / 30 // W 36" / 48 / 60 / 72



FLEX METAL TABLE SYSTEM




Height adjustable // 30-36"
 Fixed style or mobile base
 4" heavy duty locking casters optional
 Service supply modules available for easy access to h/c water, gas, electricity & data
 Under-bench storage options // drawer, cupboard or combination storage
 Adjustable reagent shelving add-ons

PRODUCT NUMBER	TABLETOP
QSKY2424	24"D X 24"W
QSKY2430	24"D X 30"W
QSKY2436	24"D X 36"W
QSKY2448	24"D X 48"W
QSKY2460	24"D X 60"W
QSKY2472	24"D X 72"W
QSKY2496	24"D X 96"W
// FOR HEIGHT ADJUSTABLE FLEX	

ALSO AVAILABLE D 30"/ 36"/ 48"


CLASSIC METAL / HEIGHT ADJUSTABLE TABLE



Painted steel frame construction
 Adjustable height legs // minimal tools required for readjustment
 Weight capacity of table // w/ levelers 600 lb. // w/ casters 300 lb.

PRODUCT NUMBER	TABLETOP
QS950H3600	27 7/8"D X 36"W
QS950H3020	27 7/8"D X 48"W
QS950H3030	27 7/8"D X 60"W
QS950H3040	27 7/8"D X 72"W
// FOR CLASSIC W/ LEVELERS	

LEVELING CASTER KIT



Set of 4 lockable casters
 Impact resistant black rubber wheels
 Mount height 3 9/16"H

PRODUCT NUMBER	HEIGHT
QS950H3200BLK	3 9/16"