



PLANNING

AXIS WORKSTATIONS



APPLICATIONS

CHEMISTRY / BIOLOGY / S.T.E.M.

FURNITURE LIST

Axis Workstations / Teacher Demonstration Desk / Wall Mounted
Drying Rack / Stainless Utility Sink / Eyewash-Safety Shower / Wood
Casework / Air Foil Fume Hoods

AXIS WORKSTATION FEATURES

8' WORK SURFACE

HEIGHT ADJUSTABLE - 30" TO 36"

PLUMBING, ELECTRIC & DATA SUPPLY CAPABLE


FEDERAL ADA GUIDELINE COMPLIANT FOR ACCESSIBILITY

19-1/2 GALLON CANTILEVERED SINK OPTIONAL

UP TO 12 ELECTRICAL OUTLETS AND DATA PORTS

VANDAL-RESISTANT UNICAST FIXTURES // CW/HW/G


EPOXY RESIN WORKTOP WITH A MARINE EDGE




SHOWN *FIXED AXIS W/ SINK*
ADDITIONAL PRODUCTS *page 32*
WOOD CASEWORK
STANDARD DEMONSTRATION DESK
TALL STORAGE CABINET
AIR FOIL FUME HOODS




AXIS WORKSTATION OPTIONS



FIXED



FIXED + SINK



MOBILE