



# PLANNING

## TEII WORK CENTERS



### APPLICATIONS

**CHEMISTRY / BIOLOGY / S.T.E.M.**

### FURNITURE LIST

TEII Work Centers / Demonstration Desk / Wall Mounted  
Drying Rack / Perimeter Sink / Eyewash-Safety Shower / Wood  
Casework / Fume Hoods

### TEII WORK CENTER FEATURES

FIRE-RETARDANT UNDERSTRUCTURE  
10-GAUGE STEEL SUPPORT PEDESTAL  
MARINE EDGE TOP TO CONTROL SPILLS  
DRAWER, CUPBOARD OR COMBINATION UNDER STORAGE  
VANDAL-RESISTANT UNICAST FIXTURES // CW/G  
RAISED TURRETS FOR TALL BEAKERS OR SUPPLIES  
ELECTRICAL OPTIONS AVAILABLE



**SHOWN** *SPLIT TOP TEII W/ FULL TROUGH*

**ADDITIONAL PRODUCTS** page 32


*CASEWORK*

*DEMONSTRATION DESK*


*TALL STORAGE CABINET*

*FUME HOODS*

TEII WORK CENTER OPTIONS



SPLIT TOP/ FULL TROUGH



FLAT TOP/ CENTER SINK