

# TEII

## WORKCENTERS

THE **TEII WORKCENTER** provides each student with ample front-counter space and elbow room. Other workers and activities are always within the student's angle of vision, making lab work not only easier but safer. Each student can work with either neighbor across from his or her corner, services such as water and electrical are easily accessible, and optional storage is located for immediate access.

### FEATURES

- + Fire-Retardant Under Structure Made of Molded Fiberglass [bolted to a 10-gauge steel support pedestal]
- + Vandal-Resistant Unicast Fixtures [CW/G]
- + Raised Turrets [for tall beakers or supplies]
- + Marine Edge Top [to control spills]
- + Optional Under Storage [drawer, cupboard or combination]
- + Electrical Options Available

